

## Human DERA Recombinant Protein

Catalog Number: PROTQ9Y315

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Overview

|              |                                                                                                                       |
|--------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name | Human DERA Recombinant Protein                                                                                        |
| Description  | Human DERA Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-318aa. Application(s): SDS-PAGE. |
| Size         | Starting from 10ug                                                                                                    |
| Tag          | His-Tag                                                                                                               |
| Form         | Liquid                                                                                                                |
| Source       | E. coli                                                                                                               |
| Formulation  | 20mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 0.1M NaCl, 1mM DTT                                             |

### Concentration

0.5mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, Q9Y315, 1-318aa; MGSSHHHHHH SSGLVPRGSH MSAHNRGTEL DLSWISIKV NHPAVLRRAE QIQARRTVKK EWQAAWLLKA VTFIDLTTLS GDDTSSNIQR LCYKAKYPIR EDLLKALNMH DKGITTAAVC VYPARVCDV KALKAAGCNI PVASVAAGFP AGQTHLKTRL EEIRLAVEDG ATEIDVVINR SLVLTGQWEA LYDEIRQFRK ACGEAHLKTI LATGELGTLT NvyKASMIAM MAGSDFIKTS TGKETVNATF PVAIVMLRAI RDFFWKTGNK IGFKPAGGIR SAKDSLAWLS LVKEELGDEW LKPELFRIGA STLLSDIERQ IYHHVTGRYA AYHDLPMs

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

### Human DERA Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.