

TELO2 (NM_016111) Human Recombinant Protein

Catalog Number: PROTQ9Y4R8

Overview

| | |
|--------------|---|
| Product Name | TELO2 (NM_016111) Human Recombinant Protein |
| Description | Recombinant protein of human TEL2, telomere maintenance 2, homolog (<i>S. cerevisiae</i>) (TELO2) |
| Size | 20 µg |
| Tag | C-Myc/DDK |
| Form | Frozen Solution in PBS Buffer |
| Source | HEK293T |
| Formulation | 25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol |

Concentration

>50 µg/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

MEPAPSEVRLAVREAIHALSSSEDGGHIFCTLES LKRYLGEMEPPALPREKEEFASAHFSPVLRCLASRLSPA WLELLPHGRLEELWASFFLEGPADQAFLVL
METIEGAAGPSFRLMKMARLLARFLREGRLAVLME AQCRQQTRPGFILLRETLGKVVALPDHLGNRLQ QENLAIEFFPQNYFRLLGEEVVRVLQAVVDS
LQGGLDSSVSFVSQVLGKACVHGRQQEILGVLV PRLAALTQGSYLHQVRCWRLVEQVPDRAMEAVL TGLVEAALGPEVLSRLLGNLVVKNKKAQFVMT
QKLLFLQSRLLTPMLQSLGHLAMDSQRRPLLQVL KELLETTWSSSAIRHTPLPQQRHVSKAVLICLA QLGEPELRDSRDELLASMMAGVKCRLDSSLPP
VRRLLGMIVAEVVSARIHPEGPLKQYEEDELSLE LLALASPQAGDGASEAGTSLVPATAEPPAETPAE IVDGGVPQAQLAGSDSDLSDDEFVYPDMS
GDRELKSSKAPAYVRDCVEALTTSEIERWEAALR ALEGLVYRSPTATREVSVELAKVLLHLEEKTCV VGFAGLRQRALVAVTVTDPAPVADYLT SQFYAL
NYSLRQRMDILDVLTAAQELSRPGCLGRTPQPGSP SPNTPCLPEAAVSQPGSAVASDWRVVVEERIRSKT QRLSKGGPRQGPAGSPSRFNSVAGHFFF P
LLQRFDRLPVTFDILLGEDQLVLGRLAHTLGALM CLAVNTTVAVAMGKALLEFVWALRFHIDAYVRQGLLS AVSSVLLSLPAARLLEDLMDELLEARSWLA
DVAEKDPDED CRTLALRALLLQLRKNRLLPPASP

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



TELO2 (NM_016111) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.