

Human UCH-L5/UCH37 Recombinant Protein

Catalog Number: PROTQ9Y5K5

Overview

Product Name	Human UCH-L5/UCH37 Recombinant Protein
Description	Human UCH-L5/UCH37 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-329aa. Application(s): Enzyme Activity, SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 5mM DTT, 30% glycerol, 200mM NaCl, 0.1mM PMSF, 2mM EDTA

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q9Y5K5, 1-329aa; MGSSHHHHHH SSSLVPRGSH MTGNAGEWCL MESDPGVFTE LIKGFGRGA QVEEWSLEP ENFEKLPVH
GLIFLKWQP GEEPAGSVVQ DSRLDTIFFA KQVINNACAT QAI SVLLNC THQDVHLGET LSEFKFSQS FDAAMKGLAL SNSDVIRQVH
NSFARQQMFE FDTKTSKKEE DAFHFVSYVP VNGRLYELDG LREGPIDLGA CNQDDWISAV RPIEKRIQK YSEGEIRFNL MAIVSDRKMI
YEQKIAELQR QLAE EEPMDT DQGNSMLSAI QSEVAKNQML IEEVQKLKR YKIENIRRKH NYLPFIMELL KTLAEHQQLI PLVEKAKEKQ
NAKKAQETK

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human UCH-L5/UCH37 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.