

STON1 (NM_006873) Human Recombinant Protein

Catalog Number: PROTQ9Y6Q2

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	STON1 (NM_006873) Human Recombinant Protein
Description	Recombinant protein of human stonin 1 (STON1)
Size	20 µg
Tag	C-Myc/DDK
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Concentration

>50 ug/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

MCSTNPGKWWTFDDDDPAVQSSQKSKNFLENQGVCRPNGLKLNLPGLREFPSGSSSTSSTPLSSPIVDFYFSPGPPSNSPLSTPTKDFPGFPGIPKAGTHV
LYPIESSSDSPLAISGGESELLPTRPTCLSHALLPSDHSCTHPTPKVGLPDEVNPQQAESLGFQSDDLQFQYFREDCAFSSPFWKDEGSDSHFTLDPPGS
KKMFSSRNKEMPIDQKSLNKCSLNYICEKLEHLQSAENQDSLRLSMHCLCAEENASSFVPHTLFRSQPKSGWSFMLRIPEKKNMMSSRQWGPFLKVL
PGGILQMYEEQGLEKPFKEIQLDPYCRLSEPKVENFSVAGKIHTVKIEHVSYTEKRKYHSKTEVVHEPDIEQMLKLGSTSYHDFLDFTTVEEELMKLPAVS
KPKKNYEEQEISLEIVDNFWGKVTKEGKFVESAVITQMYCLCFVNGNLECFLLNDLELPKRDESYEYKDKSEKKGIDILDYHFHKCVNVQFEQSRRIKFV
PLDACRFELMRFKTLYNGDNLPSLSKVVVVQGAYVELQAFVNMASLAQRSSYAGSLRSCDNIRIHFVPVPSQWIKALWTMNLQRQKSLKAKMNRACL

GSLQELESEPVIQVTVGSAKYESAYQAVVWKIDRLPDKNSSLDHPHCLSYKLELGSDQEIPSDWYPFATVQFSVPDTCASRTEVRS LGVESDVQPQKHV
QQRACYNIQVEIEKKWIKIDGEDPDKIGDCITQ

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

STON1 (NM_006873) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.