

Human DUSP10 Recombinant Protein

Catalog Number: PROTQ9Y6W6

Overview

Product Name	Human DUSP10 Recombinant Protein
Description	Human DUSP10 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 149-482aa. Application(s): SDS-PAGE, Enzyme Activity.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 2mM DTT, 50% glycerol, 200mM NaCl

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q9Y6W6, 149-482aa; I IYPN DLAKKMTKCS KSHLPSQGPV IIDCRPFMEY NKSHIQGAVH INCADKISRR RLQQGKITVL DLISCREGKD SFKRIFSKKEI IVYDENTNEP SRVMPSQPLH IVLES LKREG KEPLVLKGG L SSFKQNHENL CDNSLQLQEC REVGGGASAA SLLPQPIPT TPD IENAELT PILPFLFLGN EQDAQDLDTM QRLNIGYVIN VTTHLPLYHY EKGLFNYKRL PATDSNKQNL RQYFEEAFEF IEEAHQCGKG LLIHCQAGVS RSATIVIAYL MKHTRMTMTD AYKFVKGKRP IISPNLNFMG QLLEFEEDLN NGVTPRILTP KLMGVETVV

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human DUSP10 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.