

## Human ICOS Recombinant Protein

Catalog Number: PROTQ9Y6W8

### Overview

Product Name	Human ICOS Recombinant Protein
Description	Human ICOS Recombinant Protein expressed in Baculovirus with hlgG-His-tag. Sequence domain: 21-140aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 4ug
Tag	hlgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	20mM MES buffer (pH 5.5) containing 40% glycerol, 2mM DTT, 1mM EDTA 0.1M NaCl

### Concentration

0.25mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, Q9Y6W8, 21-140aa; EINGSAN YEMFIFHNGG VQILCKYPDI VQQFKMQLLK GGQILCDLTK TKGSGNTVSI KSLKFCHSQL SNNSVSFFLY NLDHSHANY Y FCNLSIFDPP PFKVTLTGGY LHIYESQLCC QLK

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human ICOS Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.