

Mouse recombinant Thrombospondin 4/THBS4 protein

Catalog Number: PROTQ9Z1T2

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Mouse recombinant Thrombospondin 4/THBS4 protein
Description	Mouse recombinant Thrombospondin 4/THBS4 protein expressed in HEK293 with His-Tag. Sequence domain: 27-963aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 0.1mM PMSF, 20% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Mouse, Q9Z1T2, 27-963aa; QATPQVFDLL PSSSQRLNPS ALQPVLTDPT LHEVYLISTF KLQSKSSATI FGLYSSSDNS KYFEFTVMGR LNKAILRYLK NDGKIHLVVF NNLQLADGRR HRVLLRLSNL QRGDGSVELY LDCAQADSVR NLPRAFSGLT QNPESIELRT FQRKPQDFLE ELKLVRGSL FQVASLQDCF LQQSEPLAAT STGDFNRQFL GQMTQLNQLL GEVKDLLRQQ VKETSFLRNT IAECQACGPL SFQSPTPNTL VPIAPPAPT RPTRHCDSSP CFRGVRCTDT RDGFQCGPCP DGYTGNGITC SDVDECKYHP CYPGVRCVNL APGFRCDACP VGFTGPMVQG VGINFAKTNK

QVCTDVDECQ NGACVLNSIC INTLGSYRCG PCKPGYTG DQ TRGCKTERSC RNPEQNPCSV HAQCIEERQG DVTCVCGVGW AGDGYVCGKD
VDIDSY PDEE LPCSARNCKK DNCKYVPNSG QEDADRDGIG DACDEDADGD GILNEQDNCV LTHNIDQRNS DKDIFGDACD NCRMVLNNDQ
KDTDGDGRGD ACDDDM DGDG IKNILDNCPR VPNRDQQDRD GDDVG DACDS CPDVS NPNQS DVDNDLVGDS CDTNQDSDGD
GHQDSTDNCP TVINSSQLDT DKDGIGDECD DDDDNDGIPD LVPPGPDNCR LVPNPAQEDS NNDGVGDICE ADFDQDQVID HIDVCPENAE
ITLTDFRAYQ TVVLDPEGDA QIDPNWVVLN QGMEIVQTMN SDPGLAVGYT AFNGVDFEGT FHVNTQTDDD YAGFIFGYQD SSSFYVVMWK
QTEQTYWQAT PFRAVAEPGI QLKAVKSKTG PGEHLRNSLW HTGDTSDQVR LLWKDSRNVG WKDKVSYRWF LQHRPQVG YI RVRFYEGSEL
VADSGVTIDT TMRGGRLGVF CFSQENIIWS NLKYRCNDTI PEDFQEFQTQ SFDRLDN

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Mouse recombinant Thrombospondin 4/THBS4 protein](#)

For Research Use Only. Not for use in diagnostic procedures.