

Mouse IL-10 Recombinant Protein

Catalog Number: R00021-3

Overview

Product Name	Mouse IL-10 Recombinant Protein
Description	IL-10 is an anti-inflammatory cytokine and a member of the IL-10 family of cytokines, which have indispensable functions in many infectious and inflammatory diseases. Mouse IL-10 Recombinant Protein is purified interleukin-10 produced in yeast.
Size	5 µg/vial
Tag	
Form	Lyophilized
Source	Yeast
Formulation	Lyophilized without carrier protein.

Storage

Store at -20°C for one year. Avoid repeated freeze-thaw cycles.

Purity

Ion-exchange chromatography.

Biological Activity and Protein Authenticity

In testing.

Amino Acid Sequence

SRGQYSREDN NCTHFPVGQS HMLLELR TAF SQVKTF FQTK DQLDNILLTD SLMQDFKGYL GCQALSEMIQ FYLVEVMPQA EKHGPEIKEH
LNSLGEK LKT LRMRLRRCHR FLPCENKSKA VEQVKSDFNK LQDQGVYKAM NEFDIFINCI EAYMMIKMKS (160)

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse IL-10 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.