

Mouse VEGF-A Recombinant Protein

Catalog Number: R00045-4

Overview

Product Name	Mouse VEGF-A Recombinant Protein
Description	VEGF-A is member of the family of vascular endothelial growth factor (VEGF) proteins, which stimulate vasculogenesis and angiogenesis to restore oxygen supply to tissues. VEGF-A is also a vasodilator and increases microvascular permeability and was originally referred to as vascular permeability factor (VPF). Mouse VEGF-A Recombinant Protein is purified vascular endothelial growth factor A produced in yeast.
Size	5 µg/vial
Tag	
Form	Lyophilized
Source	Yeast
Formulation	Lyophilized without carrier protein.

Storage

Store at -20°C for one year. Avoid repeated freeze-thaw cycles.

Purity

Ion-exchange chromatography.

Biological Activity and Protein Authenticity

In testing.

Amino Acid Sequence

APTTEGEQKS HEVIKFMVDVY QRSYCRPIET LVDIFQEYPD EIEYIFKPSC VPLMRCAGCC NDEALECVPT SESNITMQIM RIKPHQSQHI
GEMSFLQHSR CECRPKKDRT KPNHCEPCS ERRKHLFVQD PQTCKCCKN TDSRCKARQL ELNERTCRCD KPRR (164)

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse VEGF-A Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.