

Mouse IL-13 Recombinant Protein

Catalog Number: R00077-4

Overview

Product Name	Mouse IL-13 Recombinant Protein
Description	IL-13 is an important mediator of allergic inflammation and disease. In addition to effects on immune cells, IL-13 is implicated as a central mediator of the physiologic changes induced by allergic inflammation in many tissues. Mouse IL-13 Recombinant Protein is purified interleukin-13 produced in yeast.
Size	5 µg/vial
Tag	
Form	Lyophilized
Source	Yeast
Formulation	Lyophilized without carrier protein.

Storage

Store at -20°C for one year. Avoid repeated freeze-thaw cycles.

Purity

Ion-exchange chromatography.

Biological Activity and Protein Authenticity

In testing.

Amino Acid Sequence

GPVPRSVSLP LTLKELIEEL SNITQDQTPL CNGSMWWSVD LAAGGFCVAL DSLTNISNCN AIYRTQRILH GLCNRKAPTT VSSLPDTKIE VAHFITKLLS
YTKQLFRHGP F (111)

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse IL-13 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.