

Mouse IL-1 beta Recombinant Protein

Catalog Number: R00101-7

Overview

Product Name	Mouse IL-1 beta Recombinant Protein
Description	IL-1 beta (IL-1beta) is a member of the interleukin 1 family of cytokines. The IL-1 beta cytokine is produced by activated macrophages as a proprotein, which is proteolytically processed to its active form by caspase 1 (CASP1/ICE). This cytokine is an important mediator of the inflammatory response, and is involved in a variety of cellular activities, including cell proliferation, differentiation, and apoptosis. Mouse IL-1 beta Recombinant Protein is purified interleukin-1 beta cytokine produced in yeast.
Size	5 µg/vial
Tag	
Form	Lyophilized
Source	Yeast
Formulation	Lyophilized without carrier protein.

Storage

Store at -20°C for one year. Avoid repeated freeze-thaw cycles.

Purity

Ion-exchange chromatography.

Biological Activity and Protein Authenticity

In testing.

Amino Acid Sequence

VPIRQLHYRL RDEQQKSLVL SDPYELKALH LNGQNIQVQV IFSMSFVQGE PSNDKIPVAL GLKGKNLYLS CVMKDGTPTL QLESVDPKQY
PKKKMEKRFV FNKIEVSKV EFESAEPNW YISTSQAEHK PVFLGNNSGQ DIIDFTMESV SS (152)

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse IL-1 beta Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.