

# Human IL-6 Recombinant Protein

Catalog Number: R00102-5

## Overview

Product Name	Human IL-6 Recombinant Protein
Description	Interleukin-6 (IL-6) is an interleukin that acts as both a pro-inflammatory and anti-inflammatory cytokine. Human IL-6 Recombinant Protein is purified interleukin-6 produced in yeast.
Size	5 µg/vial
Tag	
Form	Lyophilized
Source	Yeast
Formulation	Lyophilized without carrier protein.

## Concentration

## Storage

Store at -20°C for one year. Avoid repeated freeze-thaw cycles.

## Purity

Ion-exchange chromatography.

## Biological Activity and Protein Authenticity

In testing.

## Amino Acid Sequence

VPPGEDSKDV AAPHRQPLTS SERIDKQIRY ILDGISALRK ETCNKSNMCE SSKEALAENN LNLPKMAEKD GCFQSGFNEE TCLVKIITGL LEFEVYLEYL QNRFESSEEQ ARAVQMSTKV LIQFLQKKAK NLDAITTPDP TTNASLLTKL QAQNQWLQDM TTHLILRSFK EFLQSSLRAL RQM (183)

/?/description\_after\_attributes?

## Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human IL-6 Recombinant Protein