

## Mouse IL-6 Recombinant Protein

Catalog Number: R00102-6

### Overview

Product Name	Mouse IL-6 Recombinant Protein
Description	Interleukin-6 (IL-6) is an interleukin that acts as both a pro-inflammatory and anti-inflammatory cytokine. Mouse IL-6 Recombinant Protein is purified interleukin-6 produced in yeast.
Size	5 µg/vial
Tag	
Form	Lyophilized
Source	Yeast
Formulation	Lyophilized without carrier protein.

### Concentration

### Storage

Store at -20°C for one year. Avoid repeated freeze-thaw cycles.

### Purity

Ion-exchange chromatography.

### Biological Activity and Protein Authenticity

In testing.

### Amino Acid Sequence

FPTSQVRRGDFTEDTTPNRPVYTTSQVGGILITHVLWEIVEMRKELCNGNSDCMNNDDALA  
ENNLKLPEIQRNDGCYQTGYNQEICLLKISSGLLEYHSYLEYMKNNLKDKNKKDKARVLQR  
DTETLIHIFNQEVKDLHKIVLPTPISNALLTDKLESQKEWLRTKTIQFILKSLEEFVKVT LRSTRQT (187)

🔗/description\_after\_attributes🔗

## Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse IL-6 Recombinant Protein