

Mouse IL-18 Recombinant Protein

Catalog Number: R00124-2

Overview

Product Name	Mouse IL-18 Recombinant Protein
Description	Interleukin-18 is a proinflammatory cytokine that belongs to the IL-1 superfamily. IL-18 works together with IL-12 to induce cell-mediated immunity following infection, and IL-18 is also able to induce severe inflammatory reactions, which suggests its role in certain inflammatory disorders. Mouse IL-18 Recombinant Protein is purified interleukin-18 produced in yeast.
Size	5 µg/vial
Tag	
Form	Lyophilized
Source	Yeast
Formulation	Lyophilized without carrier protein.

Storage

Store at -20°C for one year. Avoid repeated freeze-thaw cycles.

Purity

Ion-exchange chromatography.

Biological Activity and Protein Authenticity

In testing.

Amino Acid Sequence

NFGRLHCTTA VIRNINDQVL FVDKRQPVFE DMTDIDQSAS EPQTRLIIYM YKDSEVRGLA VTLSVKDSKM STLSCKNKII SFEEMDPEN
IDDIQSDLIF FQKRVPGHNK MEFESSLYEG HFLACQKEDD AFKLILKKKD ENGDKSVMFT LTNLHQ (157)

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse IL-18 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.