

Mouse CXCL10 Recombinant Protein

Catalog Number: R00278-5

Overview

Product Name	Mouse CXCL10 Recombinant Protein
Description	The ELR-negative CXC chemokine CXCL10 (IP-10) has been attributed to several roles, such as chemoattraction for monocytes/macrophages, T cells, NK cells, and dendritic cells, promotion of T cell adhesion to endothelial cells, antitumor activity, and inhibition of bone marrow colony formation and angiogenesis. Mouse CXCL10 Recombinant Protein is purified CXCL10 (IP-10) produced in yeast.
Size	5 µg/vial
Tag	
Form	Lyophilized
Source	Yeast
Formulation	Lyophilized without carrier protein.

Storage

Store at -20°C for one year. Avoid repeated freeze-thaw cycles.

Purity

Ion-exchange chromatography.

Biological Activity and Protein Authenticity

In testing.

Amino Acid Sequence

IPLARTVRCN CIHIDDGPVR MRAIGKLEII PASLSCRVE IATMKKNDE QRCLNPESKT IKNLMKAFSQ KRKRAP (77)

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse CXCL10 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.