

Mouse TNFSF11 Recombinant Protein

Catalog Number: R00363

Overview

Product Name	Mouse TNFSF11 Recombinant Protein
Description	Tumor necrosis factor ligand superfamily member 11 (TNFSF11), also known as Receptor Activator of Nuclear factor Kappa-B Ligand (RANKL), is a member of the tumor necrosis factor (TNF) ligand family. RANKL regulates bone remodeling and the development of the thymus, lymph nodes and mammary glands. Mouse TNFSF11 (RANKL) Recombinant Protein is purified TNFSF11 (RANKL) produced in yeast.
Size	5 µg/vial
Tag	
Form	Lyophilized
Source	Yeast
Formulation	Lyophilized without carrier protein.

Storage

Store at -20°C for one year. Avoid repeated freeze-thaw cycles.

Purity

Ion-exchange chromatography.

Biological Activity and Protein Authenticity

In testing.

Amino Acid Sequence

FSGAPAMMEG SWLDVAQRGK PEAQPFAHLT INAASIPSGS HKVTLSSWYH DRGWAKISNM TLSNGKLRVN QDGFYLYAN ICFRHHETSG SVPTDYQLQM VYVVKTSIKI PSSHNLMKGG STKNWSGNSE FHFYSINVGG FFKLRAGEEI SIQVSNPSLL DPDQDATYFG AFKVQDID (178)

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse TNFSF11 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.