

Mouse IL-2 Recombinant Protein

Catalog Number: R00387-6

Overview

Product Name	Mouse IL-2 Recombinant Protein
Description	Interleukin-2 (IL-2) is a cytokine produced by T-helper cells in response to antigenic or mitogenic stimulation. It is required for T-cell proliferation and other activities crucial to the regulation of the immune response. Mouse IL-2 Recombinant Protein is purified interleukin-2 produced in yeast.
Size	5 µg/vial
Tag	
Form	Lyophilized
Source	Yeast
Formulation	Lyophilized without carrier protein.

Storage

Store at -20°C for one year. Avoid repeated freeze-thaw cycles.

Purity

Ion-exchange chromatography.

Biological Activity and Protein Authenticity

In testing.

Amino Acid Sequence

APTSSSTSSP TSSSTAEAAQQ QQQQQQQQQQ HLEQLLMDLQ ELLSRMENYR NLKLPRLMTF KFYLPKQATE LKDLQCLEDE LGPLRHVLDL
TQSKSFQLED AENFISNIRV TVVKLKGSDN TFECQFDDES ATVVDLRRW IAFCSIIST SPQ (153)

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse IL-2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.