

## Mouse CCL3 Recombinant Protein

Catalog Number: R00405-2

### Overview

Product Name	Mouse CCL3 Recombinant Protein
Description	Chemokine ligand 3 (CCL3), also known as Macrophage inflammatory protein-1 alpha (MIP-1 alpha) is a small cytokine belonging to the CC chemokine family. CCL3/MIP-1 alpha is a chemoattractant for several different leukocytes, with varying degrees of potency (monocytes > T-cells > neutrophils > eosinophils).
Size	5 µg/vial
Tag	
Form	Lyophilized
Source	Yeast
Formulation	Lyophilized without carrier protein.

### Storage

Store at -20°C for one year. Avoid repeated freeze-thaw cycles.

### Purity

Ion-exchange chromatography.

### Biological Activity and Protein Authenticity

In testing.

### Amino Acid Sequence

APFGADTPTA CCFSYVARQL SRKIVADYFE TSSQCSKPGV IFQTKKGRQV CANPTEDWVQ EYITDLELNA (70)

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Mouse CCL3 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.