

Human IL-17A Recombinant Protein

Catalog Number: R00421-3

Overview

Product Name	Human IL-17A Recombinant Protein
Description	IL-17A is a member of the IL-17 family of cytokines, whose members are involved in numerous immune regulatory functions. IL-17 induces the production of many other cytokines, chemokines, and prostaglandins. Human IL-17A Recombinant Protein is purified interleukin-17A produced in yeast.
Size	5 µg/vial
Tag	
Form	Lyophilized
Source	Yeast
Formulation	Lyophilized without carrier protein.

Storage

Store at -20°C for one year. Avoid repeated freeze-thaw cycles.

Purity

Ion-exchange chromatography.

Biological Activity and Protein Authenticity

In testing.

Amino Acid Sequence

GITIPRNPGC PNSDKNFPR TVMVNLNIHN RNTNTNPKRS SDYYNRSTSP WNLHRNEDPE RYPSVIWEAK CRHLGCINAD GNVDYHMNSV
PIQQEILVLR REPPHCPNSF RLEKILVSVG CTCVTPIVHH VA (132)

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human IL-17A Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.