

## Mouse IL-17A Recombinant Protein

Catalog Number: R00421-4

### Overview

Product Name	Mouse IL-17A Recombinant Protein
Description	IL-17A is a member of the IL-17 family of cytokines, whose members are involved in numerous immune regulatory functions. IL-17 induces the production of many other cytokines, chemokines, and prostaglandins. Mouse IL-17A Recombinant Protein is purified interleukin-17A produced in yeast.
Size	5 µg/vial
Tag	
Form	Lyophilized
Source	Yeast
Formulation	Lyophilized without carrier protein.

### Storage

Store at -20°C for one year. Avoid repeated freeze-thaw cycles.

### Purity

Ion-exchange chromatography.

### Biological Activity and Protein Authenticity

In testing.

### Amino Acid Sequence

AAIIPQSSAC PNTEAKDFLQ NVKVNLIKVFN SLGAKVSSRR PSDYLNRRSTS PWTLHRNEDP DRYPSVIWEA QCRHQRCVNA EGKLDHHMNS  
VLIQEQEILVL KREPESCPFT FRVEKMLVGV GCTCVASIVR QAA (133)

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Mouse IL-17A Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.