

## Mouse CCL4 Recombinant Protein

Catalog Number: R00703-4

### Overview

Product Name	Mouse CCL4 Recombinant Protein
Description	The small chemotactic cytokine CCL4 (MIP-1 beta) is a chemoattractant for natural killer cells, monocytes and a variety of other immune cells. CCL4 is involved in several inflammatory and autoimmune diseases including viral infection such as HIV-1/AIDS. Mouse CCL4 Recombinant Protein is purified chemokine ligand 4 (CCL4, MIP-1 beta) produced in yeast.
Size	5 µg/vial
Tag	
Form	Lyophilized
Source	Yeast
Formulation	Lyophilized without carrier protein.

### Storage

Store at -20°C for one year. Avoid repeated freeze-thaw cycles.

### Purity

Ion-exchange chromatography.

### Biological Activity and Protein Authenticity

In testing.

### Amino Acid Sequence

APMGSDPPTS CCFSYTSRQL HRSFVMDYYE TSSLCSKPAV VFLTKRGRQI CANPSEPWVT EYMSDLELN (69)

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Mouse CCL4 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.