

Rat CCL4 Recombinant Protein

Catalog Number: R00703-6

Overview

Product Name	Rat CCL4 Recombinant Protein
Description	The small chemotactic cytokine CCL4 (MIP-1 beta) is a chemoattractant for natural killer cells, monocytes and a variety of other immune cells. CCL4 is involved in several inflammatory and autoimmune diseases including viral infection such as HIV-1/AIDS. Rat CCL4 Recombinant Protein is purified chemokine ligand 4 (CCL4, MIP-1 beta) produced in yeast.
Size	5 µg/vial
Tag	
Form	Lyophilized
Source	Yeast
Formulation	Lyophilized without carrier protein.

Storage

Store at -20°C for one year. Avoid repeated freeze-thaw cycles.

Purity

Ion-exchange chromatography.

Biological Activity and Protein Authenticity

In testing.

Amino Acid Sequence

APIGSDPPTS CCFSYTSRKI HRNFVMDYYE TSSLCSQPAV VFLTKKGRQI CADPSEPWVN EYVNDLELN (69)

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Rat CCL4 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.