

Mouse GM-CSF Recombinant Protein

Catalog Number: R00911-5

Overview

Product Name	Mouse GM-CSF Recombinant Protein
Description	Granulocyte-macrophage colony-stimulating factor (GM-CSF) is a protein secreted by macrophages, T cells, mast cells, NK cells, endothelial cells and fibroblasts. GM-CSF stimulates stem cells to produce granulocytes (neutrophils, eosinophils, and basophils) and monocytes. Mouse GM-CSF Recombinant Protein is purified GM-CSF produced in yeast.
Size	5 µg/vial
Tag	
Form	Lyophilized
Source	Yeast
Formulation	Lyophilized without carrier protein.

Storage

Store at -20°C for one year. Avoid repeated freeze-thaw cycles.

Purity

Ion-exchange chromatography.

Biological Activity and Protein Authenticity

In testing.

Amino Acid Sequence

APTRSPITVT RPWKHVEAIK EALNLLDDMP VTLNEEVEVV SNEFSFKKLT CVQTRLKIFE QGLRGNFTKL KGALNMTASY YQTYCPPTPE
TDCETQVTY ADFIDSLKTF LTDIPFECKK PGQK (124)

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse GM-CSF Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.