

Mouse BAFF Recombinant Protein

Catalog Number: R01257-3

Overview

Product Name	Mouse BAFF Recombinant Protein
Description	BAFF (B-cell Activation Factor), also known as Tumor Necrosis Factor Superfamily member 13B (TNFSF13B), is expressed in B cell lineage cells and acts as a potent B cell activator. It has been also shown to play an important role in the proliferation and differentiation of B cells. Mouse BAFF Recombinant Protein is purified b-cell activation factor produced in yeast.
Size	5 µg/vial
Tag	
Form	Lyophilized
Source	Yeast
Formulation	Lyophilized without carrier protein.

Storage

Store at -20°C for one year. Avoid repeated freeze-thaw cycles.

Purity

Ion-exchange chromatography.

Biological Activity and Protein Authenticity

In testing.

Amino Acid Sequence

AFQGPETEIQ DVDLSAPPAP CLPGCRHSQH DDNGMNLRLNI IQDCLQLIAD SDTPTIRKGT YTFVPWLLSF KRGNAL EEKE NKIVVRQTGY
FFIYSQVLYT DPIFAMGHVI QRKKVHVFGD ELSLVTLFRC IQNMPKTLPN NSCYSAGIAR LEEGDEIQLA IPRENAQISR NGDDTFFGAL KLL (183)

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse BAFF Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.