

Mouse APRIL Recombinant Protein

Catalog Number: R02417-2

Overview

Product Name	Mouse APRIL Recombinant Protein
Description	A proliferation-inducing ligand (APRIL), also known as TNFSF13, is a member of the tumor necrosis factor (TNF) ligand superfamily. APRIL/TNFSF13 has been shown to play a role in protecting cells from undergoing apoptosis and promoting B cell development. Mouse APRIL Recombinant Protein is purified APRIL (TNFSF13) produced in yeast.
Size	5 µg/vial
Tag	
Form	Lyophilized
Source	Yeast
Formulation	Lyophilized without carrier protein.

Storage

Store at -20°C for one year. Avoid repeated freeze-thaw cycles.

Purity

Ion-exchange chromatography.

Biological Activity and Protein Authenticity

In testing.

Amino Acid Sequence

AVLTQKHKKK HSVLHLVPVN ITSKADSDVT EVMWQPVLRG GRGLEAQGDI VRVWDTGIYLLYSQVLFHDV TFTMGQVVSRE GQGRRETLF
RCIRSMPSDP DRAYNSCYSA GVFHLHQGDI ITVKIPRANA KLSLSPHGTF LGFVKL (146)

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse APRIL Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.