

## Human IL-19 Recombinant Protein

Catalog Number: R04434

### Overview

Product Name	Human IL-19 Recombinant Protein
Description	Interleukin-19 (IL-19) is a cytokine that belongs to the IL-10 family. It is up-regulated in monocytes following stimulation with granulocyte-macrophage colony-stimulating factor (GM-CSF) or lipopolysaccharide. Human IL-19 Recombinant Protein is purified interleukin-19 produced in yeast.
Size	5 µg/vial
Tag	
Form	Lyophilized
Source	Yeast
Formulation	Lyophilized without carrier protein.

### Storage

Store at -20°C for one year. Avoid repeated freeze-thaw cycles.

### Purity

Ion-exchange chromatography.

### Biological Activity and Protein Authenticity

In testing.

### Amino Acid Sequence

LRRLISTDM HHIEESFQEI KRAIQAKDTF PNVTLSTLE TLQIIKPLDV CCVTKNLLAF YVDRVFKDHQ EPNPKILRKI SSIANSFLYM QKTLRQCQEQ  
RQCHCRQEAT NATRVIHDNY DQLEVHAAAI KSLGELDVFL AWINKNHEVM FSA (153)

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human IL-19 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.