

Human VEGF-D Recombinant Protein

Catalog Number: R04809

Overview

Product Name	Human VEGF-D Recombinant Protein
Description	VEGF-D is member of the family of vascular endothelial growth factor (VEGF) proteins, which stimulate vasculogenesis and angiogenesis to restore oxygen supply to tissues. Human VEGF-D Recombinant Protein is purified vascular endothelial growth factor D produced in yeast.
Size	5 µg/vial
Tag	
Form	Lyophilized
Source	Yeast
Formulation	Lyophilized without carrier protein.

Storage

Store at -20°C for one year. Avoid repeated freeze-thaw cycles.

Purity

Ion-exchange chromatography.

Biological Activity and Protein Authenticity

In testing.

Amino Acid Sequence

FAATFYDIET LKVIDEEWQR TQCSPRETCV EVASELGKST NTFFKPPCVN VFRCGGCCNE ESLICMNTST SYISKQLFEI SVPLTSVPEL
VPVKVANHTG CKCLPTAPRH PYSIIRR (117)

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human VEGF-D Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.