

Mouse IL-34 Recombinant Protein

Catalog Number: R06903

Overview

Product Name	Mouse IL-34 Recombinant Protein
Description	Interleukin-34 (IL-34) is a homodimeric cytokine recognized by colony-stimulating factor 1 receptor (CSF-1R). IL-34 is involved in the survival, proliferation, and differentiation of mononuclear phagocytes. Mouse IL-34 Recombinant Protein is purified interleukin-15 produced in yeast.
Size	5 µg/vial
Tag	
Form	Lyophilized
Source	Yeast
Formulation	Lyophilized without carrier protein.

Storage

Store at -20°C for one year. Avoid repeated freeze-thaw cycles.

Purity

Ion-exchange chromatography.

Biological Activity and Protein Authenticity

In testing.

Amino Acid Sequence

NENLEIWTLT QDKECDLTGY LRGKLYKYNR LQYMKHYFPI NYRIAVPYEG VLRVANITRL QKAHVSEREL RYLWVLSLN ATESVMDVLL
EGHPSWKYLQ EVQTLLENVQ RSLMDVEIGP HVEAVLSLLS TPGLSLKLVK PKALLDNCFR VMELLYCSCC KQSPILKWQD CELPRLPHS
PGSLMQCTAT NVYPLSRQTP TSLPGSPSS HGSLP (215)

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse IL-34 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.