

SARS-CoV-2 (COVID-19) NSP5A Protein, His Tag

Catalog Number: RCOV01

Overview

Product Name	SARS-CoV-2 (COVID-19) NSP5A Protein, His Tag
Description	SARS-CoV-2 (COVID-19) NSP5A Protein, His Tag is produced in E. coli and has a theoretical molecular weight of 34.6KD. This product is for research use only. Product is under validation for additional applications and indications. If you're interested, please contact support@bosterbio.com.
Size	100ug/vial
Tag	
Form	Lyophilized
Source	E. coli
Formulation	Lyophilized from sterile 20mM PB, 150mM NaCl, PH 7.3-7.4, 10% glycerol and 4% trehalose

Storage

The product is shipped at ambient temperature. Upon receipt, store it immediately at -20°C for 6 months under sterile conditions. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Purity

> 90%. Purification measurement method by SDS-PAGE quantitative densitometry by Coomassie® Blue Staining.

Method of purification: Nickel column affinity purification

Amino Acid Sequence

6×His tag at C-terminal

Accession #: YP_009725301

SGFRKMAFPSGKVEGCMVQVTCGTTTLNGLWLDDVVYCPRHVICTSEDMLNPYEDLLIRKSNHNFLVQAGNVQLRVIGHSMQNCVLKLVDTANP
KTPKYKfVRIQPGQTfSVLACYNGSPSGVYQCAMRPNFTIKGSFLNGSCGSVGFNIDYDCVSFCYMHMELPTGVHAGTDLEGNFYGPfVDRQTAQA
AGTDTTITVNVLAWLYAAVINGDRWFLNRFtTTLNDFNLVAMKYNIEPLTQDHVDILGPLSAQTGIAVLDMCASLKELLQNGMNGRTILGSALLEDEFT
PFDVVRQCSGVTFQ

Background

Coronaviruses (CoV) are a family of large and enveloped positive-sense single-stranded RNA viruses that are classified into four genera, the alpha, beta, gamma, and delta coronaviruses. While gamma and delta coronaviruses mainly infect birds, alpha and beta coronaviruses are known to infect mammals and cause human respiratory illnesses such as the common cold, pneumonia, and severe diseases like SARS, MERS, and COVID-19.

Coronavirus nucleoproteins localize to the cytoplasm and the nucleolus, a subnuclear structure, in both virus-infected primary cells and in cells

transfected with plasmids that express N protein. Coronavirus N protein is required for coronavirus RNA synthesis, and has RNA chaperone activity that may be involved in template switch. Nucleocapsid protein is the most abundant protein of coronavirus. During virion assembly, the N protein binds to viral RNA and leads to the formation of the helical nucleocapsid. Nucleocapsid protein is a highly immunogenic phosphoprotein also implicated in viral genome replication and in modulating cell signaling pathways. Because of the conservation of N protein sequence and its strong immunogenicity, the N protein of coronavirus is chosen as a diagnostic tool.

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



SARS-CoV-2 (COVID-19) NSP5A Protein, His Tag

For Research Use Only. Not for use in diagnostic procedures.