

Anti-Fibroblast growth factor 1 FGF1 Antibody Picoband®

Catalog Number: RP1005

About FGF1

Fibroblast growth factor 1 (acidic), also known as FGF1/ECGF/HBGF1, is a human gene which is mapped to 5q31. Human FGF1 shares 96% amino acid sequence homology with both rat and mouse. The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities, and are involved in a variety of biological processes, including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth and invasion. This protein functions as a modifier of endothelial cell migration and proliferation, as well as an angiogenic factor. It acts as a mitogen for a variety of mesoderm- and neuroectoderm-derived cells in vitro, thus is thought to be involved in organogenesis. Additionally, Acidic fibroblast growth factor is derived from beta-endothelial cell growth factor (ECGFB) by posttranslational processing. Alpha-ECGF is also derived from ECGFB in the same manner.

Overview

Product Name	Anti-Fibroblast growth factor 1 FGF1 Antibody Picoband®
Reactive Species	Human
Description	Boster Bio Anti-Fibroblast growth factor 1 FGF1 Antibody catalog # RP1005. Tested in ELISA, IHC, WB applications. This antibody reacts with Human. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	ELISA, IHC, WB
Clonality	Polyclonal
Formulation	Each vial contains 0.9mg NaCl, 0.2mg Na ₂ HPO ₄ , 0.05mg Na ₃ . Carrier free (No BSA) form available in stock. If you want this antibody carrier free please specify "Carrier Free" or "No BSA" in your order note.
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P05230

Technical Details


Immunogen	E. coli-derived human FGF1 recombinant protein (Position: F16-D155).
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot, and HRP Conjugated anti-Rabbit IgG Super Vision Assay Kit (SV0002-1) for IHC(P).
Cross Reactivity	No cross-reactivity with other proteins

Isotype	Rabbit IgG
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Immunohistochemistry (Paraffin-embedded Section), 0.5-1ug/ml, Human ELISA, 0.1-0.5ug/ml, - Western blot, 0.1-0.5ug/ml, Human

Anti-Fibroblast growth factor 1 FGF1 Antibody Picoband® (RP1005) Images



Anti-human FGF1 antibody, RP1005, IHC(P)IHC(P): Human Mammary Cancer Tissue



Anti-human FGF1 antibody, RP1005, Western blotting
Lane 1: Recombinant Human FGF1 Protein 10ng
Lane 2: Recombinant Human FGF1 Protein 5ng
Lane 3: Recombinant Human FGF1 Protein 2

2 Publications Citing This Product

1. PubMed ID: 25667662, Zhang L, Wang H, Wang T, Jiang N, Yu P, Chong Y, Fu F. Exp Ther Med. 2015 Mar;9(3):972-976. Epub 2014 Dec 24. Ferulic Acid Ameliorates Nerve Injury Induced By Cerebral Ischemia In Rats.
2. PubMed ID: 21461910, Zhou Y, Ren L, Zhu J, Yan S, Wang H, Song N, Li L, Ouyang H, Pang D. Mol Cell Biochem. 2011 Aug;354(1-2):39-46. Doi: 10.1007/S11010-011-0803-8. Epub 2011 Apr 3. Construction Of A Recombinant Human Fgf1 Expression Vector For Mammary Gland-Specific ...

Visit bosterbio.com/anti-human-fgf1-antibody-rp1005-boster.html to see all 2 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Fibroblast growth factor 1 FGF1 Antibody

For Research Use Only. Not for use in diagnostic procedures.