

Anti-Interleukin-10 IL10 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: RP1015-Fluoro647

About IL10

Interleukin-10 (IL-10 or IL10), also known as human cytokine synthesis inhibitory factor (CSIF), is an anti-inflammatory cytokine. In humans IL-10 is encoded by the IL10 gene. It is capable of inhibiting synthesis of pro-inflammatory cytokines like IFN-gamma, IL-2, IL-3, TNFalpha and GM-CSF made by cells such as macrophages and regulatory T-cells. IL-10 also displays potent abilities to suppress the antigen presentation capacity of antigen presenting cells. Kim et al. (1992) showed that the mouse IL 10 gene contains 5 exons and spans about 5.2 kb of genomic DNA. Eskdale et al. (1997) mapped the IL10 gene to the junction between 1q31 and 1q32.

Overview

Product Name	Anti-Interleukin-10 IL10 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Mouse
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P18893

Technical Details

Immunogen	E. coli-derived mouse IL-10 recombinant protein (Position: S19-S178).
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

19 Publications Citing This Product

1. PubMed ID: 10.3389/fmicb.2020.622354, Prevention and Alleviation of Dextran Sulfate Sodium Salt-Induced Inflammatory Bowel Disease in Mice With Bacillus subtilis-Fermented Milk via Inhibition of the Inflammatory Responses and Regulation of the Intestinal Flora
2. PubMed ID: 10.1016/j.ajpath.2021.09.006, Chronic intestinal inflammation suppresses brain activity by inducing neuroinflammation in mice
3. PubMed ID: 10.1007/s00436-020-07025-3, The role of regulatory B cells in Echinococcus granulosus -infected mice

Visit bosterbio.com/anti-mouse-il-10-antibody-rp1015-boster.html to see all 19 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Interleukin-10 IL10 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.