

## Anti-Interleukin-10 IL10 Antibody Picoband® HRP Conjugated

Catalog Number: RP1015-HRP

### About IL10

Interleukin-10 (IL-10 or IL10), also known as human cytokine synthesis inhibitory factor (CSIF), is an anti-inflammatory cytokine. In humans IL-10 is encoded by the IL10 gene. It is capable of inhibiting synthesis of pro-inflammatory cytokines like IFN-gamma, IL-2, IL-3, TNFalpha and GM-CSF made by cells such as macrophages and regulatory T-cells. IL-10 also displays potent abilities to suppress the antigen presentation capacity of antigen presenting cells. Kim et al. (1992) showed that the mouse IL 10 gene contains 5 exons and spans about 5.2 kb of genomic DNA. Eskdale et al. (1997) mapped the IL10 gene to the junction between 1q31 and 1q32.

### Overview

Product Name	Anti-Interleukin-10 IL10 Antibody Picoband® HRP Conjugated
Reactive Species	Mouse
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na2HPO4.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P18893

### Technical Details

Immunogen	E. coli-derived mouse IL-10 recombinant protein (Position: S19-S178).
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

### 19 Publications Citing This Product

1. PubMed ID: 10.3389/fmicb.2020.622354, Prevention and Alleviation of Dextran Sulfate Sodium Salt-Induced Inflammatory Bowel Disease in Mice With Bacillus subtilis-Fermented Milk via Inhibition of the Inflammatory Responses and Regulation of the Intestinal Flora
2. PubMed ID: 10.1016/j.ajpath.2021.09.006, Chronic intestinal inflammation suppresses brain activity by inducing neuroinflammation in mice
3. PubMed ID: 10.1007/s00436-020-07025-3, The role of regulatory B cells in Echinococcus granulosus -infected mice

Visit [bosterbio.com/anti-mouse-il-10-antibody-rp1015-boster.html](https://bosterbio.com/anti-mouse-il-10-antibody-rp1015-boster.html) to see all 19 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Interleukin-10 IL10 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.