

## Anti-IL12 p40/IL12B Antibody Picoband® Biotin Conjugated

Catalog Number: RP1028-Biotin

### About IL12B

Subunit beta of interleukin 12, also known as interleukin 12B, or IL12B, is a subunit of human interleukin 12. It is mapped to 5q33.3. This cytokine is expressed by activated macrophages that serve as an essential inducer of Th1 cells development. This cytokine has been found to be important for sustaining a sufficient number of memory/effector Th1 cells to mediate long-term protection to an intracellular pathogen. Overexpression of this gene was observed in the central nervous system of patients with multiple sclerosis (MS), suggesting a role of this cytokine in the pathogenesis of the disease.

### Overview

Product Name	Anti-IL12 p40/IL12B Antibody Picoband® Biotin Conjugated
Reactive Species	Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9R278

### Technical Details

Immunogen	E.coli-derived rat IL-12 (p40) recombinant protein (Position: M23-S321).
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

1. PubMed ID: 19538477, Anti-arthritis activity of cationic materials

Visit [bosterbio.com/anti-il-12-p40-antibody-rp1028-boster.html](http://bosterbio.com/anti-il-12-p40-antibody-rp1028-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL12 p40/IL12B Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.