

## Anti-Claudin 2/CLDN2 Antibody Picoband® APC Conjugated

Catalog Number: RP1042-APC

### About CLDN2

Claudin-2 is a protein that in humans is encoded by the CLDN2 gene. By genomic sequence analysis, this gene is mapped to chromosome Xq22.3-q23. This gene product belongs to the claudin protein family whose members have been identified as major integral membrane proteins localized exclusively at tight junctions. Claudins are expressed in an organ-specific manner and regulate tissue-specific physiologic properties of tight junctions. By expression in a human intestinal epithelial cell line, it was determined that the intestine-specific homeodomain proteins CDX1 and CDX2 activated a reporter plasmid driven by the CLDN2 promoter.

### Overview

Product Name	Anti-Claudin 2/CLDN2 Antibody Picoband® APC Conjugated
Reactive Species	Human, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P57739

### Technical Details

Immunogen	E.coli-derived human Claudin 2 recombinant protein (Position: A38-V230). Human Claudin 2 shares 92% amino acid (aa) sequence identity with mouse Claudin 2.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

---

1. PubMed ID: 23600383, Cen J, Liu L, Li Ms, He L, Wang Lj, Liu Yq, Liu M, Ji Bs. J Pharm Pharmacol. 2013 May;65(5):665-72. Doi: 10.1111/Jphp.12033. Epub 2013 Feb 26. Alteration In P-Glycoprotein At The Blood-Brain Barrier In The Early Period Of Mcao In Rats.

Visit [bosterbio.com/anti-claudin-2-antibody-rp1042-boster.html](http://bosterbio.com/anti-claudin-2-antibody-rp1042-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Claudin 2/CLDN2 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.