

Anti-Transcription factor MafA MAFA Antibody Picoband®

Catalog Number: RP1066

About MAFA

MAFA is a transcription factor that binds RIPE3b, a conserved enhancer element that regulates pancreatic beta cell-specific expression of the insulin gene. MAFA plays a key role in coordinating and controlling the level of insulin gene expression in islet beta cells. It is located in 8q24.3. The intracellular domain of MAFA contains a single immunoreceptor tyrosine-based inhibitory motif (ITIM), which classifies MAFA as a member of an inhibitory receptor superfamily. The inhibitory capacity of MAFA is best defined in mast cells, where MAFA keeps in check the antigen-induced (i.e. Fc epsilon RI-induced) secretion of inflammatory mediators.

Overview


Product Name	Anti-Transcription factor MafA MAFA Antibody Picoband®
Reactive Species	Human, Rat
Description	Boster Bio Anti-Transcription factor MafA MAFA Antibody catalog # RP1066. Tested in WB applications. This antibody reacts with Human, Rat. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	WB
Clonality	Polyclonal
Formulation	Each vial contains antibody formulated with stabilizing components, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄ , and 0.05 mg NaN ₃ . *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8NHW3

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human MAFA, identical to the related mouse sequence.
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot.

Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.1-0.5ug/ml, Human, Rat

Anti-Transcription factor MafA MAFA Antibody Picoband® (RP1066) Images



Western blot analysis of MAFA using anti-MAFA antibody (RP1066). Electrophoresis was performed on a 5-20% SDS-PAGE gel at 70V (Stacking gel) / 90V (Resolving gel) for 2-3 hours. The sample well of each lane was loaded with 50ug of sample under reducing conditions. Lane 1: human HEK293 whole cell lysates, Lane 2: human K562 whole cell lysates, Lane 3: human 22RV1 whole cell lysates, Lane 4: rat PC-12 whole cell lysates. After Electrophoresis, proteins were transferred to a Nitrocellulose membrane at 150mA for 50-90 minutes. Blocked the membrane with 5% Non-fat Milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-MAFA antigen affinity purified polyclonal antibody (Catalog # RP1066) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an Enhanced Chemiluminescent detection (ECL) kit (Catalog # EK1002) with Tanon 5200 system. A specific band was detected for MAFA at approximately 37KD. The expected band size for MAFA is at 37KD.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Transcription factor MafA MAFA Antibody

For Research Use Only. Not for use in diagnostic procedures.