

## Anti-Transcription factor MafA MAFA Antibody Picoband® Biotin Conjugated

Catalog Number: RP1066-Biotin

### About MAFA

MAFA is a transcription factor that binds RIPE3b, a conserved enhancer element that regulates pancreatic beta cell-specific expression of the insulin gene. MAFA plays a key role in coordinating and controlling the level of insulin gene expression in islet beta cells. It is located in 8q24.3. The intracellular domain of MAFA contains a single immunoreceptor tyrosine-based inhibitory motif (ITIM), which classifies MAFA as a member of an inhibitory receptor superfamily. The inhibitory capacity of MAFA is best defined in mast cells, where MAFA keeps in check the antigen-induced (i.e. Fc epsilon RI-induced) secretion of inflammatory mediators.

### Overview

Product Name	Anti-Transcription factor MafA MAFA Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8NHW3

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human MAFA, identical to the related mouse sequence.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Transcription factor MafA MAFA Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.