

## Anti-Actin/ACTA1 Antibody Picoband® APC Conjugated

Catalog Number: RP1070-APC

### About ACTA1

Actin, a highly conserved protein, is a major component of both the cytoskeletal and contractile structures in the cell types. It varies in amount, being related to the type of differentiation and to the functional state of cells and tissues. The actins exhibit over 90% sequence homology, but each isoform has a unique NH<sub>2</sub>-terminal sequence. The isoforms are comprised of three alpha-actin, one beta-actin, two gamma-actin. Because the amino acid sequence of the C-terminal is the same for almost all actins, this antibody has been raised using a synthetic peptide corresponding to the C-terminal 11 residues.

### Overview

Product Name	Anti-Actin/ACTA1 Antibody Picoband® APC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P68133

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human Actin, identical to the related mouse and rat sequences.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## 65 Publications Citing This Product

1. PubMed ID: 32390601, Cheng J,Luan J,Chen P,Kuang X,Jiang P,Zhang R,Chen S,Cheng F,Gou X.Immunosuppressive receptor LILRB1 acts as a potential regulator in hepatocellular carcinoma by integrating with SHP1.Cancer Biomark.2020; 28(3):309-319.doi:10.3233/CBM-190940.PMID:32390601

2. PubMed ID: 33221387, Sandulescu T, Weniger J, Philippou S, Mücke T, Naumova EA, Arnold WH. Immunohistochemical evidence of striated muscle cells within midfacial superficial musculoaponeurotic system. Ann Anat. 2020 Nov 19:151647. doi: 10.1016/j.aanat.2020.151647. Epub ahead

3. PubMed ID: 15237432, Inhibition of human La protein by RNA interference downregulates hepatitis B virus mRNA in 2.2.15 cells

Visit [bosterbio.com/anti-actin-antibody-rp1070-boster.html](https://bosterbio.com/anti-actin-antibody-rp1070-boster.html) to see all 65 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Actin/ACTA1 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.