

## Anti-JNK2/MAPK9 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: RP1077-Fluoro594

### About MAPK9

JNK2 is also known as MAPK9. The protein encoded by this gene is a member of the MAP kinase family. MAP kinases act as an integration point for multiple biochemical signals, and are involved in a wide variety of cellular processes such as proliferation, differentiation, transcription regulation and development. This kinase targets specific transcription factors, and thus mediates immediate-early gene expression in response to various cell stimuli. It is most closely related to MAPK8, both of which are involved in UV radiation induced apoptosis, thought to be related to the cytochrome c-mediated cell death pathway. Also, this gene and MAPK8 are also known as c-Jun N-terminal kinases. This kinase blocks the ubiquitination of tumor suppressor p53, and thus it increases the stability of p53 in nonstressed cells. Studies of this gene's mouse counterpart suggest a key role in T-cell differentiation. Several alternatively spliced transcript variants encoding distinct isoforms have been reported.

### Overview

Product Name	Anti-JNK2/MAPK9 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P45984

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human JNK2, identical to the related mouse and rat sequences.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

## 2 Publications Citing This Product

1. PubMed ID: 29393346, Wang P, Mao Z, Pan Q, Lu R, Huang X, Shang X, Zhang R, You H. Int J Mol Med. 2018 Apr;41(4):2117-2127. doi: 10.3892/ijmm.2018.3410. Epub 2018 Jan 22. Histone deacetylase-4 and histone deacetylase-8 regulate interleukin-1beta-induced cartilage catabol...

2. PubMed ID: 27293989, Recombinant Newcastle disease virus (rL-RVG) triggers autophagy and apoptosis in gastric carcinoma cells by inducing ER stress

Visit [bosterbio.com/anti-jnk2-antibody-rp1077-boster.html](http://bosterbio.com/anti-jnk2-antibody-rp1077-boster.html) to see all 2 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-JNK2/MAPK9 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.