

Anti-Osteopontin/SPP1 Antibody Picoband® HRP Conjugated

Catalog Number: RP1080-HRP

About SPP1

Osteopontin (OPN), also known as secreted phosphoprotein 1 (SPP1), is a protein that in humans is encoded by the SPP1 gene (secreted phosphoprotein 1). The protein encoded by this gene is involved in the attachment of osteoclasts to the mineralized bone matrix. And the encoded protein is secreted and binds hydroxyapatite with high affinity. The osteoclast vitronectin receptor is found in the cell membrane and may be involved in the binding to this protein. Also, this protein is a cytokine that upregulates expression of interferon-gamma and interleukin-12. Several transcript variants encoding different isoforms have been found for this gene.

Overview

Product Name	Anti-Osteopontin/SPP1 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P10451

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human Osteopontin, different from the related mouse sequence by eight amino acids, and from the related rat sequence by seven amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

10 Publications Citing This Product

1. PubMed ID: 33650080, Li X,Hu B,Wang L,Xia Q,Ni X. P2X7 receptor-mediated phenotype switching of pulmonary artery smooth muscle cells in hypoxia. Mol Biol Rep.2021 Mar 1.doi:10.1007/s11033-021-06222-2. Epub ahead of print.PMID:33650080.
2. PubMed ID: -, Caiyun Mu,Ye He,Yan Hu,Menghuan Li,Maowen Chen,Rong Wang,Yang Xiang, Zhong Luo,Kaiyong Cai,Construction of chemokine substance P-embedded biomimetic multilayer onto bioactive magnesium silicate-titanium implant for bone regeneration,Applied Materials Today
3. PubMed ID: -, Chaowei Hu,Kun Zuo,Kuibao Li,Yuanfeng Gao,Mulei Chen,Roumu Hu,Ye Liu,Hongjie Chi,Hongjiang Wang,Yanwen Qin,Xiaoyan Liu, Jiuchang Zhong,Jun Cai,Xinchun Yang,Jing Li, "p38/JNK Is Required for the Proliferation and Phenotype Changes of Vascular Smooth Muscle

Visit bosterbio.com/anti-osteopontin-antibody-rp1080-boster.html to see all 10 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Osteopontin/SPP1 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.