

Anti-Protein Wnt-7a Wnt7a Antibody Picoband®

Catalog Number: RP1088

About WNT7A

This gene is a member of the WNT gene family, which consists of structurally related genes that encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is involved in the development of the anterior-posterior axis in the female reproductive tract, and also plays a critical role in uterine smooth muscle patterning and maintenance of adult uterine function. Mutations in this gene are associated with Fuhrmann and Al-Awadi / Raas - Rothschild / Schinzel phocomelia syndromes.

Overview


Product Name	Anti-Protein Wnt-7a Wnt7a Antibody Picoband®
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Protein Wnt-7a Wnt7a Antibody catalog # RP1088. Tested in WB applications. This antibody reacts with Human, Mouse, Rat. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	WB
Clonality	Polyclonal
Formulation	Each vial contains 4 mg Trehalose, 0.9 mg NaCl and 0.2 mg Na ₂ HPO ₄ .
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O00755

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human Wnt7a, identical to the related mouse sequence.
Recommended Detection Systems	Boster recommends ECL Plus Western Blotting Substrate (Catalog # AR1196-200) for Western blot.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Lyophilized

Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.1-0.5ug/ml, Human, Mouse, Rat

Anti-Protein Wnt-7a Wnt7a Antibody Picoband® (RP1088) Images



Western blot analysis of Wnt7a using anti-Wnt7a antibody (PB9198). Electrophoresis was performed on a 10% SDS-PAGE gel at 80V (Stacking gel) / 120V (Resolving gel) for 2 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: human hepatocellular carcinoma paracancerous tissue (HCCP) lysates, Lane 2: rat lung tissue lysates, Lane 3: mouse lung tissue lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-Wnt7a antigen affinity purified polyclonal antibody (PB9198) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody (Catalog # BA1054) at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an ECL Plus Western Blotting Substrate (Catalog # AR1196-200) with Tanon 5200 system. A specific band was detected for Wnt7a at approximately 39 kDa. The expected band size for Wnt7a is at 39 kDa.

1 Publications Citing This Product

1. PubMed ID: 33359712, An R,Wang X,Yang L,Zhang J,Wang N,Xu F,Hou Y,Zhang H,Zhang L. Polystyrene microplastics cause granulosa cells apoptosis and fibrosis in ovary through oxidative stress in rats.Toxicology.2020 Dec 24;449:152665. doi:10.1016/j.tox.2020.152665.Epub ahead of p

Visit bosterbio.com/anti-wnt7a-antibody-rp1088-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Protein Wnt-7a Wnt7a Antibody

For Research Use Only. Not for use in diagnostic procedures.