

## Anti-Macrophage metalloelastase MMP12 Antibody Picoband® APC Conjugated

Catalog Number: RP1089-APC

### About Mmp12

Matrix metalloproteinase-12 (MMP12), also known as MME or ME, is an enzyme that in humans is encoded by the MMP12 gene. The gene is part of a cluster of MMP genes which localize to chromosome 11q22.2. It is thought that the protein encoded by this gene is cleaved at both ends to yield the active enzyme, but this processing has not been fully described. The enzyme degrades soluble and insoluble elastin. It may play a role in aneurysm formation and studies in mice suggest a role in the development of emphysema. This gene may be involved in tissue injury and remodeling.

### Overview

Product Name	Anti-Macrophage metalloelastase MMP12 Antibody Picoband® APC Conjugated
Reactive Species	Mouse
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P34960

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of mouse MMP12, different from the related human sequence by thirteen amino acids, and from the related rat sequence by three amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

1. PubMed ID: 25489966, Vaz M, Rajasekaran S, Potteti Hr, Reddy Sp. Am J Respir Cell Mol Biol. 2015 Jul;53(1):125-34. Doi: 10.1165/Rcmb.2014-0118Oc. Myeloid-Specific Fos-Related Antigen-1 Regulates Cigarette Smoke-Induced Lung Inflammation, Not Emphysema, In Mice.

Visit [bosterbio.com/anti-mmp12-antibody-rp1089-boster.html](http://bosterbio.com/anti-mmp12-antibody-rp1089-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Macrophage metalloelastase MMP12 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.