

Anti-CD153/TNFSF8 Antibody

Catalog Number: RP1103

About TNFSF8

CD153 is a ligand for CD30. The protein encoded by this gene is a cytokine that belongs to the tumor necrosis factor (TNF) ligand family. This cytokine is a ligand for TNFRSF8/CD30, which is a cell surface antigen and a marker for Hodgkin lymphoma and related hematologic malignancies. The engagement of this cytokine expressed on B cell surface plays an inhibitory role in modulating Ig class switch. This cytokine was shown to enhance cell proliferation of some lymphoma cell lines, while to induce cell death and reduce cell proliferation of other lymphoma cell lines. The pleiotropic biologic activities of this cytokine on different CD30+ lymphoma cell lines may play a pathophysiologic role in Hodgkin's and some non-Hodgkin's lymphomas. Two transcript variants encoding different isoforms have been found for this gene.

Overview


Product Name	Anti-CD153/TNFSF8 Antibody
Reactive Species	Human
Description	Boster Bio Anti-CD153/TNFSF8 Antibody catalog # RP1103. Tested in Flow Cytometry, IHC applications. This antibody reacts with Human.
Application	Flow Cytometry, IHC
Clonality	Polyclonal
Formulation	Each vial contains antibody formulated with stabilizing components, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄ , and 0.05 mg NaN ₃ . *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P32971

Technical Details


Immunogen	E.coli-derived human CD153 recombinant protein (Position: Q63-D234). Human CD153 shares 75% amino acid (aa) sequence identity with mouse CD153.
Recommended Detection Systems	Boster recommends HRP Conjugated anti-Rabbit IgG Super Vision Assay Kit (SV0002-1) for IHC(P).
Cross Reactivity	No cross-reactivity with other proteins.

Isotype	Rabbit IgG
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Immunohistochemistry (Paraffin-embedded Section), 0.5-1ug/ml Flow Cytometry, 1-3ug/1x10 ⁶ cells

Anti-CD153/TNFSF8 Antibody (RP1103) Images



CD153 was detected in paraffin-embedded sections of human tonsil tissues using rabbit anti-CD153 Antigen Affinity purified polyclonal antibody (Catalog # RP1103) at 1 ug/mL. The immunohistochemical section was developed using SABC method (Catalog # SA1022).



2. Flow Cytometry analysis of U937 cells using anti-CD153 antibody (RP1103). Overlay histogram showing U937 cells stained with RP1103 (Blue line). The cells were blocked with 10% normal goat serum. And then incubated with rabbit anti-CD153 Antibody (RP1103, 1ug/1x10⁶ cells) for 30 min at 20°C. DyLight488 conjugated goat anti-rabbit IgG (BA1127, 5-10ug/1x10⁶ cells) was used as secondary antibody for 30 minutes at 20°C. Isotype control antibody (Green line) was rabbit IgG (1ug/1x10⁶) used under the same conditions. Unlabelled sample (Red line) was also used as a control.

1 Publications Citing This Product

1. PubMed ID: 29984229, CD30 Is Highly Expressed in Chronic Obstructive Pulmonary Disease and Induces the Pulmonary Vascular Remodeling

Visit bosterbio.com/anti-cd153-antibody-rp1103-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD153/TNFSF8 Antibody

For Research Use Only. Not for use in diagnostic procedures.