

Anti-SIAH Interacting Protein/CACYBP Antibody Picoband® Fluoro550 Conjugated

Catalog Number: RP1104-Fluoro550

About CACYBP

Calcyclin-binding protein is a protein that in humans is encoded by the CACYBP gene. And this gene is mapped to 1q24-q25. The protein encoded by this gene is a calyculin binding protein. It may be involved in calcium-dependent ubiquitination and subsequent proteosomal degradation of target proteins. In addition, it probably serves as a molecular bridge in ubiquitin E3 complexes and participates in the ubiquitin-mediated degradation of beta-catenin. Two alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Overview

Product Name	Anti-SIAH Interacting Protein/CACYBP Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, IP, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9HB71

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human CACYBP, different from the related mouse sequence by five amino acids, and from the related rat sequence by six amino acids.
Cross Reactivity	No cross-reactivity with other proteins
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SIAH Interacting Protein/CACYBP Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.