

Cy3 Conjugated AffiniPure Rabbit Anti-Goat IgG (H+L)

Basic Information

Catalog Number	BA1034
Size	0.5 mL
Concentration	1 mg/mL
Ig Type	IgG
Immunogen	Goat IgG (whole molecular)
Host	Rabbit
Purity	The antibody was purified from antiserum by immunoaffinity chromatography.
Conjugation	Cy3
Fluorophore	Amax=554nm; Emax=568nm
Form	Concentrated, Liquid
Tested Application	IF, Flow Cytometry
Contents	0.5 mg of affinity purified rabbit anti-goat IgG (Purify>99%), F/P=3-9/1, 0.01 M PBS, 5 mg/mL BSA, 50% glycerol.
Specificity	This Cy3 conjugated antibody is specific for goat IgG and shows no cross-reactivity with human/rat/mouse/rabbit IgG.

Suggested Working Concentration

Immunofluorescence: 2-10 µg/mL

Flow Cytometry: 2-10 µg/mL

Optimal working dilutions must be determined by end users.

Application Notes

Please centrifuge before use.

It is concentrated solution, should be diluted by neutral PBS or TBS before use.

The working reagent (WR) after diluting should be used up within a day, otherwise, please give up the surplus.

Storage

At -20° C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.