

# Fluoro594 Conjugated AffiniPure Rabbit IgG (H+L)

## Basic Information

<b>Catalog Number</b>	BA1045-Fluoro594
<b>Size</b>	1 mL
<b>Concentration</b>	1 mg/mL
<b>Ig Type</b>	Rabbit IgG
<b>Purity</b>	The antibody was purified from antiserum by immunoaffinity chromatography.
<b>Conjugation</b>	Fluoro594
<b>Fluorophore</b>	Amax=593nm; Emax=618nm
<b>Form</b>	Concentrated, Liquid
<b>Tested Application</b>	Flow Cytometry
<b>Contents</b>	1 mg/ml of Fluoro594 conjugated specific antibody, 0.01M PBS, 5mg/ml BSA, 50% glycerol.
<b>Specificity</b>	This Fluoro594 conjugated antibody is specific for rabbit IgG.

## Suggested Working Concentration

Flow Cytometry: 2-10 µg/mL

**Optimal working dilutions must be determined by end users.**

## Application Notes

Please centrifuge before use and dilute with neutral PBS or TBS.

The diluted solution should be used on the same day. Any unused solution should be discarded.

## Storage

At -20° C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.