

PE Conjugated AffiniPure Rabbit IgG (H+L)

Basic Information

Catalog Number	BA1045-PE
Size	100 µg
Concentration	0.5 mg/mL
Ig Type	Rabbit IgG
Purity	The antibody was purified from antiserum by immunoaffinity chromatography.
Conjugation	PE
Fluorophore	Excitation Wavelength: 566nm Emission Wavelength: 574nm
Form	Concentrated, Liquid
Tested Application	Flow Cytometry
Contents	0.5 mg/ml of PE conjugated specific antibody, 0.01M PBS, 5mg/ml BSA, 50% glycerol.
Specificity	This PE conjugated antibody is specific for rabbit IgG.

Suggested Working Concentration

Flow Cytometry: 1-3 µg/1x10⁶ cells

Optimal working dilutions must be determined by end users.

Application Notes

Please centrifuge before use and dilute with neutral PBS or TBS.

Use the diluted solution immediately. Discard any unused portion.

Storage

At -20° C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.