

HRP Conjugated AffiniPure Donkey Anti-Rabbit IgG(H+L)

Basic Information

Catalog Number	BA1061
Size	100 µL
Concentration	1 mg/mL
Ig Type	IgG
Immunogen	Rabbit IgG (whole molecular)
Host	Donkey
Purity	The antibody was purified from antiserum by immunoaffinity chromatography.
Conjugation	Peroxidase
Form	Concentrated, Liquid
Tested Application	WB, ELISA
Contents	0.1 mg of HRP conjugated specific antibody, 0.01 M PBS, 50% glycerol.
Specificity	This HRP conjugated antibody is specific for Rabbit IgG and shows no cross-reactivity with Bovine, Chicken, Goat, Guinea Pig, Syrian Hamster, Horse, Human, Mouse, Rat, and Sheep Serum Proteins.

Suggested Working Concentration

Western Blotting: 0.005-1µg/ml (ECL detection)

Direct ELISA: 0.01-0.2µg/ml (TMB detection)

Optimal working dilutions must be determined by end users.

Application Notes

Dilute with neutral PBS or TBS.

Use of sodium azide as a preservative will substantially inhibit the enzyme activity of horseradish peroxidase.

Storage

At -20° C for one year from date of receipt. Avoid repeated freezing and thawing.