

DyLight 550 Conjugated AffiniPure Donkey Anti-Mouse IgG (H+L)

Basic Information

| | |
|---------------------------|---|
| Catalog Number | BA1149 |
| Size | 100 µL |
| Concentration | 1 mg/mL |
| Ig Type | IgG |
| Immunogen | Mouse IgG (whole molecular) |
| Host | Donkey |
| Purity | The antibody was purified from antiserum by immunoaffinity chromatography. |
| Conjugation | DyLight 550 |
| Fluorophore | Amax=562nm; Emax=576nm |
| Form | Concentrated, Liquid |
| Tested Application | IF, Flow Cytometry |
| Contents | 1 mg/mL of DyLight 550 conjugated specific antibody (purity is above 99%), 0.01 M PBS, 5 mg/ml BSA, 50% glycerol. |
| Specificity | This DyLight 550 conjugated antibody is specific for mouse IgG and shows no cross-reactivity with human/ bovine/rabbit IgG. |

Suggested Working Concentration

Immunofluorescence: 2-10 µg/mL

Flow Cytometry: 2-10 µg/mL

Optimal working dilutions must be determined by end users.

Application Notes

Please centrifuge before use and dilute with neutral PBS or TBS.

The diluted solution should be used on the same day. Any unused solution should be discarded.

Storage

At -20° C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.