





Product Name:	Closed System of Centrifuge Tube
Product Application:	<p>1. The pre-assembled Closed System of Centrifuge Tube is convenient for customers in biomedical/biopharmaceutical/vaccine production industries to aseptically sample and transfer the tissue and cell culture/ vaccine solutions, which can effectively save the assembling and sampling time and lower the cost.</p> <p>2. Each is double-packed and equipped with a seal cap which can be used to substitute the pipes in case of subsequent centrifugation and analysis.</p>
Product Materials:	PP (body)+ Silicone tube
product Specifications:	<p>Temperature range: stored at room temperature</p> <p>Shelf life: 3 Years after date of production (ensure package is in good)</p>
Sterilization:	<p>Yes.</p> <p>Sterilized by E-beam, Sterility Assurance Level: SAL=10⁻⁶. The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137</p>
Pyrogens:	Non-Pyrogenic
RNase/DNase:	DNase/RNase free
BSE/TSE:	These products are deemed animal free.
Performance Testing:	<p>Each manufacturing lot is sampled and tested in accordance with standard operating procedures.</p> <p>Appearance inspection: qualified Sealing test: qualified</p> <p>Packaging inspection: qualified</p>
Features	<p>1. The pre-assembled Closed System is ready to use at hand. The closed liquid transfer avoids the contamination risk of open transfer and eliminates the necessity of subsequent cleaning validation.</p> <p>2. The pipeline material meets the requirements of industrial pharmaceutical enterprises, and is connected to male taper at the connector, which prevents accidental disconnection and offers convenience.</p> <p>3. The double aseptic packaging suits GMP production.</p> <p>4. Pipelines can be customized according to demand.</p>

Product Range

Cat. No.	Product name	Picture	Product Description	/Pack	/Case
C60622-EFB050A	50mL centrifuge tubes with preassembled closed systems		<p>1. Liquid Transfer Tube</p>	1	4
C60922-EFB050A	250mL centrifuge tubes with preassembled closed systems		<p>1/8"(3.2mm) ID, 1/4"(6.4mm) OD , 50cm in length.</p> <p>2. 24mm membrane diameter; 0.2µm pore size of membrane.</p> <p>3.The end of Liquid Transfer</p>	1	4
C61022-EFB050A	500mL centrifuge tubes with preassembled closed systems		<p>Tube is connected to the male taper and Luer cap.</p>	1	4

Technical Drawing



Product picture

