



Product Name:	PCR Plate
Product Application:	For PCR machine whose fluorescent signals is read at the top, Nest white PCR plates can effectively prevent the cross interference between different holes. In addition, Nest white PCR plate can read low fluorescent signals with high stability. Nest transparent PCR plates apply to the PCR machines whose fluorescent signals is read at the bottom.
Product Materials:	Imported special PP
product Specifications:	Temperature range: stored at room temperature Shelf life: 5 years after date of production (ensure package is in good)
Sterilization:	Yes (some of all catalog numbers). Sterilized by E-beam, Sterility Assurance Level: SAL=10 ⁻⁶ . The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137
Pyrogens:	Non-Pyrogenic
RNase/DNase:	DNase/RNase free
BSE/TSE:	These products are deemed animal free.
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified Vacuum Seal test: qualified Freeze testing: qualified Centrifugation test: qualified Evaporation rate test: qualified Autoclave test : qualified Packaging inspection: qualified
Features	<ol style="list-style-type: none"> 1. Made from USP grade VI polypropylene. 1. Flat surface, thick and solid, not easy to deform. 2. The black printing mark is designed for easy identification. 3. Elevated edge of the hole can better prevent cross contamination. 4. Compatible with NEST pressure sensitive film, self-adhesive film and hot sealing film 5. Autoclavable. 6. PCR Plates with No Skirt Product Features: <ol style="list-style-type: none"> a) The low profile design minimized the condensate contamination and evaporation of reaction solution. 7. DNase/RNase free, no human DNA and PCR inhibitor. 8. PCR Plate, No Skirt, Elevated Wells Product Features: <ol style="list-style-type: none"> a) Elevated wells plate is designed to reduce the risk of cross contamination b) Easy to open 9. PCR Plates with Skirt Product Features: <ol style="list-style-type: none"> a) Stronger construction and more space for labeling.

Product Range



Cat. No.	Product Description	sterile	/Pack	/Case
402101	0.1 mL 96 Well PCR Plate, No Skirt, Clear, H12 Notch	No	25	100
402401	0.1 mL 96 Well PCR Plate, Semi Skirt, Clear, A1 Notch, Compatible with ABI Machine	No	25	100
402501	0.1 mL 96 Well PCR Plate, Full Skirt, Clear, H1 Notch	No	25	100
402511	0.1 mL 96 Well PCR Plate, Full Skirt, White, H1 Notch	No	25	100
402001	0.2 mL 96 Well PCR Plate, No Skirt, Clear, H1 Notch	No	25	100
402201	0.2 mL 96 Well PCR Plate, No Skirt, Elevated Wells, Clear, H12 Notch	No	25	100
402301	0.2 mL 96 Well PCR Plate, Semi Skirt, Clear, A12 Notch	No	25	100
402601	0.2 mL 96 Well PCR Plate, Semi Skirt, Clear, A12 Notch, Compatible with ABI Machine	No	25	100
402111	0.1 mL 96 Well PCR Plate, No Skirt, White, H12 Notch	No	25	100
402411	0.1 mL 96 Well PCR Plate, Semi Skirt, White, A1 Notch	No	25	100
402011	0.2 mL 96 Well PCR Plate, No Skirt, White, H1 Notch	No	25	100
402712	100 µl 96 Well PCR Plates, Semi Skirt, Roche, H12 Notch, White	No	10	50
402812	100 µl 96 Well PCR Plates, No Skirt, Roche, A12/H12 Notch, White	No	10	50
409013	40 µl 384 Well PCR Plates, Semi Skirt, Compatible with Roche Machine, White Frame, A24 Notch, White	No	10	50
409033 Sterilized	40 µl 384 Well PCR Plates, Semi Skirt, Compatible with Roche Machine, White Frame, A24 Notch, White. Sterilized	Yes	10	50

Technical Drawing



402001(Clear) / 402011(White)





402301(Clear)





402501(Clear)





402101(Clear) / 402111(White)





402401(Clear)/402411(White)





402201(Clear)




402601(Clear)



402812



409013



402712

